

947ml Buffered Isotonic Saline Eyewash Bottle

If your eyes become contaminated with particles or chemicals, your first response can mean the difference between being able to see or permanent eye damage.



947ml eyewash bottle
Part No. EWB947

Wall mount and bottle
holder
Part No. WMB947

Opto-Pharm Eye Wash from Spill Station is a sterile buffered isotonic solution specifically formulated to enhance eye comfort by irrigating, flushing and cleaning without interfering with the natural functions.

Opto-Pharm Eye Wash is formulated to provide optimum first response eye and wound irrigation.

The functional wall station ensures your eyewash is always nearby, accessible and easily seen.

The exclusive eyewash diffuser cup fits comfortably on the eye and ensures an even flow of flushing solution.

Simply twist the lid and start treatment.

- State of the art Blow-Fill-Seal equipment ensures safety of product.
- In house QA & QC laboratories.
- Full 4 year shelf life on Eye Wash range. (from manufacturing date) Stability Study Data available.
- Certificate of Analysis is issued with each shipment, batch specific

**Opto-Pharm Eye Wash is
manufactured to meet ANSI Z358.1-2004 (as personal eye wash only).**

